

three upstream oligonucleotides diagrammed in Figure 32 (*i.e.*, one of SEQ ID NOS:44-46), in a total volume of 5  $\mu$ l of distilled water. The reactions were overlaid with a drop of ChillOut™ evaporation barrier (MJ Research) and warmed to 62°C. The cleavage reactions were started by the addition of 5  $\mu$ l of an enzyme mixture to each tube, and the reactions were incubated at 62°C for 30 min. The reactions shown in lanes 1-3 of Figure 34 received Mixture 1; reactions 4-6 received Mixture 2; reactions 7-9 received Mixture 3 and reactions 10-12 received Mixture 4.

After 30 minutes at 62°C, the reactions were stopped by the addition of 8  $\mu$ l of 95% formamide with 20 mM EDTA and 0.05% marker dyes. Samples were heated to 75°C for 2 minutes immediately before electrophoresis through a 20% acrylamide gel (19:1 cross-linked), with 7 M urea, in a buffer of 45 mM Tris-Borate, pH 8.3, 1.4 mM EDTA.

Following electrophoresis, the products of the reactions were visualized by the use of an Hitachi FMBIO fluorescence imager, the output of which is seen in Figure 34. The reaction products shown in lanes 1, 4, 7 and 10 of Figure 34 were from reactions which contained SEQ ID NO:44 as the upstream oligonucleotide (see Figure 32a). The reaction products shown in lanes 2, 5, 8 and 11 of Figure 34 were from reactions which contained SEQ ID NO:45 as the upstream oligonucleotide (see Figure 32b). The reaction products shown in lanes 3, 6, 9 and 12 of Figure 34 were from reactions which contained SEQ ID NO:46, the invader oligonucleotide, as the upstream oligonucleotide (see Figure 32c).

Examination of the  $Mn^{2+}$  based reactions using either Cleavase® A/G nuclease or DNAPtaq as the cleavage agent (lanes 1 through 3 and 4 through 6, respectively) shows that both enzymes have active exonuclease function in these buffer conditions. The use of a 3' label on the probe oligonucleotide allows the products of the nibbling activity to remain labeled, and therefore visible in this assay. The ladders seen in lanes 1, 2, 4 and 5 confirm that the probe hybridize to the target DNA as intended. These lanes also show that the location of the non-invasive oligonucleotides have little effect on the products generated. The uniform ladder created by these digests would be difficult to distinguish from a ladder caused by a contaminating nuclease, as one

might find in a clinical specimen. In contrast, the products displayed in lanes 3 and 6, where an invader oligonucleotide was provided to direct the cleavage, show a very distinctive shift, so that the primary cleavage product is smaller than those seen in the non-invasive cleavage. This product is then subject to further nibbling in these conditions, as indicated by the shorter products in these lanes. These invader-directed cleavage products would be easily distinguished from a background of non-specific degradation of the probe oligonucleotide.

When  $Mg^{2+}$  is used as the divalent cation the results are even more distinctive. In lanes 7, 8, 10 and 11 of Figure 34, where the upstream oligonucleotides were not invasive, minimal nibbling is observed. The products in the DNAP $Taq$  reactions show some accumulation of probe that has been shortened on the 5' end by one or two nucleotides consistent with previous examination of the action of this enzyme on nicked substrates (Longley *et al.*, *supra*). When the upstream oligonucleotide is invasive, however, the appearance of the distinctively shifted probe band is seen. These data clearly indicated that it is the invasive 3' portion of the upstream oligonucleotide that is responsible for fixing the site of cleavage of the downstream probe.

Thus, the above results demonstrate that it is the presence of the free or initially non-annealed nucleotides at the 3' end of the invader oligonucleotide which mediate the shift in the cleavage site, not just the presence of an oligonucleotide annealed upstream of the probe. Nucleic acid detection assays which employ the use of an invader oligonucleotide are termed "invader-directed cleavage" assays.

#### EXAMPLE 14

Invader-Directed Cleavage Recognizes Single And Double Stranded  
Target Molecules In A Background Of Non-Target DNA Molecules

For a nucleic acid detection method to be broadly useful, it must be able to detect a specific target in a sample that may contain large amounts of other DNA, *e.g.*, bacterial or human chromosomal DNA. The ability of the invader directed cleavage

assay to recognize and cleave either single- or double-stranded target molecules in the presence of large amounts of non-target DNA was examined. In these experiments a model target nucleic acid, M13, in either single or double stranded form (single-stranded M13mp18 is available from Life Technologies, Inc and double-stranded M13mp19 is available from New England Biolabs), was combined with human genomic DNA (Novagen, Madison, WI) and then utilized in invader-directed cleavage reactions. Before the start of the cleavage reaction, the DNAs were heated to 95°C for 15 minutes to completely denature the samples, as is standard practice in assays, such as polymerase chain reaction or enzymatic DNA sequencing, which involve solution hybridization of oligonucleotides to double-stranded target molecules.

For each of the reactions shown in lanes 2-5 of Figure 35, the target DNA (25 fmole of the ss DNA or 1 pmole of the ds DNA) was combined with 50 pmole of the invader oligonucleotide (SEQ ID NO:46); for the reaction shown in lane 1 the target DNA was omitted. Reactions 1, 3 and 5 also contained 470 ng of human genomic DNA. These mixtures were brought to a volume of 10 µl with distilled water, overlaid with a drop of ChillOut™ evaporation barrier (MJ Research), and brought to 95°C for 15 minutes. After this incubation period, and still at 95°C, each tube received 10 µl of a mixture comprising 2.25 µl of Cleavase® A/G nuclease extract (prepared as described in Example 2) and 5 pmole of the probe oligonucleotide (SEQ ID NO:43), in 20 mM MOPS, pH 7.5 with 0.1 % each of Tween 20 and NP-40, 4 mM MnCl<sub>2</sub> and 100 mM KCl. The reactions were brought to 62°C for 15 minutes and stopped by the addition of 12 µl of 95% formamide with 20 mM EDTA and 0.05% marker dyes. Samples were heated to 75°C for 2 minutes immediately before electrophoresis through a 20% acrylamide gel (19:1 cross-linked), with 7 M urea, in a buffer of 45 mM Tris-Borate, pH 8.3, 1.4 mM EDTA. The products of the reactions were visualized by the use of an Hitachi FMBIO fluorescence imager. The results are displayed in Figure 35.

In Figure 35, lane 1 contains the products of the reaction containing the probe (SEQ ID NO:43), the invader oligonucleotide (SEQ ID NO:46) and human genomic DNA. Examination of lane 1 shows that the probe and invader oligonucleotides are